


Medica Display for Medical Device

# MEDDP-415

15" Slim Medical Display



## Features

- 15 " 1024 x 768 Resolution Display
- P-cap touch
- Stylish And Elegant Design

## Specifications

### DISPLAY

Size	15"
Resolution	1024 x 768
Max. Colors	16.2M
Brightness (Typ.)	300 nits
Viewing Angle (CR=10)(Typ.)	80°(H)/80°(V)
Contrast Ratio (Typ.)	700:1
Touch Screen	P-cap Touch, Protective Glass

### SYSTEM

OS Support	Windows®XP, Windows® 7, Windows® 10
------------	-------------------------------------

### I/O

Audio	Audio-in x 1
Serial Port	USB x 1 for touch controller
Signal	VGA x 1, HDMI x 1
Speaker	1 W speaker x 2 on backside

### MECHANICAL AND ENVIRONMENTAL

Power Consumption	Full Loading:24.37Watts
VESA	75 x 75mm & 100 x 100mm
Power Requirement	DC12V / 5A power input
Operating Temperature	5°C ~ 40°C(41°F ~ 104°F)
Storage Temperature	-10°C ~ 60°C(14°F ~ 140°F)
Storage Humidity	10%-90%@35°C, non-condensing
Dimension	352.5x 276.1 x 45.0 mm
Packing Size	402.5x 326.1 x 80.0 mm
Gross Weight	3.7 kg (8.16lb)
Net Weight	2.6kg (5.73 lb)
Certifications	CE: EN 60601-1-2:2015(V4.0) EN 60601-1:2006 / A1:2013(V3.1) FCC: Part 18 Class B, Part 15 Class B

## Ordering Information

- MEDDP-415HNN-A1-1010  
15" LCD Monitor.300 nits.VGAX1.HDMIx1.DC 12V Adapter
- MEDDP-415HPN-A1-1010  
15" LCD Monitor.300 nits, P-CAP touch.VGAX1.HDMIx1.DC 12V Adapter