


Medica Display for Medical Device

MEDDP-615

15.6" Slim Medical Display



Features

- 15.6" 1920 x 1080 Resolution Display
- Projected Capacitive Multi-touch
- VGA input and HDMI input
- Stylish and elegant design

Specifications

DISPLAY

| | |
|-----------------------------|------------------------------------|
| Size | 15.6" |
| Resolution | Full HD 1920 x 1080 |
| Max. Colors | 16.7M |
| Brightness (Typ.) | 300 nits |
| Viewing Angle (CR=10)(Typ.) | 85°(H)/85°(V) |
| Contrast Ratio (Typ.) | 700:1 |
| Touch Screen | Projective capacitive touch screen |

SYSTEM

| | |
|------------|--------------------------------------|
| OS Support | Windows® XP, Windows® 7, Windows® 10 |
|------------|--------------------------------------|

I/O

| | |
|-------------|------------------------------|
| Audio | Audio-in x 1 |
| Serial Port | USB x 1 for touch controller |
| Signal | VGA x 1, HDMI x 1 |
| Speaker | 1 W speaker x 2 on back side |

MECHANICAL AND ENVIRONMENTAL

| | |
|-----------------------|--|
| Power Consumption | Full Loading :24.37Watts |
| VESA | 75 x 75mm & 100 x 100mm |
| Power Requirement | DC12V / 5A power input |
| Operating Temperature | 5°C ~ 40°C(41°F ~ 104°F) |
| Storage Temperature | -10°C ~ 60°C(14°F ~ 140°F) |
| Storage Humidity | 10%~90%@35°C, non-condensing |
| Dimension | 387.8 x 232.9 x 38.4 mm |
| Packing Size | 425.3 x 269.2 x 45.6 mm |
| Gross Weight | 6.5 kg (14.33lb) |
| Net Weight | 2.8 kg (6.2lb) |
| Certifications | CE: EN 60601-1-2:2015(V4.0), EN 60950-1 EN 60601-1:2006/A1:2013/A12:2014 (V3.1) EN 62366:2007 + A1:2014 (V1.1) FCC: Part 18 Class B, Part 15 Class B UL: ANSI/AAMI ES60601-1:2012 (V3.1) cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1) |

Ordering Information

- MEDDP-615HPN-A1-1010
15.6" LCD Monitor.300 nits.P-CAP touch,VGAx1.HDMIx1.DC 12V Adapter