



## Mobile Application Computer

**Venus-123**

11.6" Fanless 6th Generation Core i5 Mobile Computer

**Features**

- Intel® 6th Generation Core i5-6300U 3.0GHz Dual Core
- 11.6" 250-nit WXGA LCD
- Supports Non-ECC DDR4 SODIMM up to 16GB
- Capacitive Multi-Touch Screen / 5-Wire Resistive Touch Screen
- Easy Cleaning Flat Screen
- 5 Programmable Smart Function Keys
- LAN/COM Medical Isolation
- Trusted Platform Module 2.0
- Remote Management: Intel Active Management Technology
- +12V DC Input
- Dual 32.7Wh Backup Batteries

**Specifications****Main Specifications**

Processor	Intel® 6th Generation Core i5-6300U 3.0GHz Dual Core
System Memory	Non-ECC DDR4 SODIMM up to 16GB
Storage Disk Drive	M.2 Solid State Disk Drive x1
OS Support	Windows 7, Window 8.1, Window 10, Linux
Wireless Communication	802.11ac/a/b/g/n, Bluetooth V4.1
Function Key	Power On/Off, Speaker Volume Up/Down, LCD Brightness Up/Down
Power Requirement	+12V DC, Dual 32.7Wh Backup Batteries
Speaker	2W x2
Security	Trusted Platform Module 2.0, RFID Reader(optional)

**Display**

Display Size	11.6" LCD
Max. Resolution	1366 x 768
Luminance(cd/m2)(TYP)	250 nits
Max. Colors	16.7 M
Viewing Angle	90°(H)/50°(V)
Contrast Ratio	500:1
Back Light Life Time	15,000 hours
Touch Screen	5-Wire Resistive / Capacitive Multi-Touch

**I/O**

USB	USB 2.0 x2, USB 3.0 x2
Serial Port	Isolated RS-232 x1
Ethernet	Isolated Gigabit LAN x 2
Vedio Out	HDMI x 1

**Mechanical and Environmental**

Power Consumption	Full loading : 29 watts
Operating Temperature	0°C ~ 35°C(32°F ~ 95°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Mounting	VESA 75/100 mm
Degree of Protection	IP54
Net Weight	1.9Kg
Gross Weight	TBD
Package Dimension (W x D x H)	536 x 230 x382 mm
Certifications	CE: EN 60601-1-2:2015(V4.0), EN 60601-1:2006/A1:2013(V3.1) FCC: Part 18 Class B UL: ANSI/AAMI ES60601-1:2012 (V3.1) cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1)

ALL PRODUCT, PRODUCT SPECIFICATIONS AND DATA ARE SUBJECT TO CHANGE WITHOUT NOTICE TO IMPROVE RELIABILITY, FUNCTION OR DESIGN OR OTHERWISE.