


Medical PC & Monitor for Digital OR

ACCEL-JS500

Nvidia Jetson AGX Xavier platform for Medical AI Application



Power button



Features

- Nvidia Jetson AGX Xavier Platform for Medical AI usage.
- 512-core Volta GPU with Tensor Cores
- Support optional 3G or 12G SDI input for image capture
- Support Dual 4K display output (HDMI and Display port)
- Support optional internal speaker
- Medical Certification
- Palm size dimension with powerful AI computing

Specifications

Main Specifications

AI Engine	Nvidia Jetson AGX Xavier
CPU	8-Core ARM v8.2 64bit CPU, 8MB L2 + 4MB L3
System Memory	32GB 256-Bit LPDDR4
Graphics	512 Core Volta GPU with Tensor Cores
Storage	32GB eMMC Onboard
	Micro SD card Slot M.2 M Key 2280
OS support	Linux with Jetpack 4.4.1
Security	TPM 2.0
Speaker	2W x Speaker x 2(Optional)

I/O


USB	USB Type C x 2 for USB 3.1
	USB Type A x 1 for USB 3.0
	USB Type A x 1 for USB 2.0
Ethernet	Gigabit LAN x 1
Video Out	HDMI 2.0 x 1, Display Port 1.4 x 1
	3G SDI x 1 or 12G SDI x 1
Video input (Optional)	3G SDI(Max resolution up to 1920x1080p@60fps)
	12G SDI(Max resolution up to 4096x2160@60fps)
Audio	Mic-in x 1, Line out x 1
COM	RS232 x 1
Function Port	Reset Button x 1, Recovery Button x 1,
DC-in	12V DC Jack
Wireless Communication	802.11a/b/g/n/ac.w/BT5.0 (Optional)
LED Indicator	Power LED x 1

Environment and Mechanical

Power	Medical adapter 12V 84W
Operating temperature	0°C ~ 35°C(32°F ~ 95°F)
Storage temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Dimension	184(W) x 145(L) x 66mm(H)
Package Size	290(W) x 225(L) x 198mm(H)
Gross Weight	TBD
Net Weight	TBD
Certifications	CE: EN 60601-1-2:2015(V4.0),
	EN 60601-1:2006/A1:2013/A12:2014 (V3.1)
	FCC: Part 18 Class B, Part 15 Class B,
	UL: ANSI/AAMI ES60601-1:2012 (V3.1)
	cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1)

Nvidia Jetson AGX Xavier platform for Medical AI Application

Dimension / Unit: mm



Ordering Information

- » ACCEL-JS500-N1-A1-0010
Medical grade AI box PC, Nvidia Jetson AGX Xavier,
Fan. adapter with internal speaker
- » ACCEL-JS500-N1-A2-0010
Medical grade AI box PC, Nvidia Jetson AGX Xavier, Fan.
3G SDI. adapter, with internal speaker
- » ACCEL-JS500-N1-A3-0010
Medical grade AI box PC, Nvidia Jetson AGX Xavier,
Fan. 12G SDI. adapter, with internal speaker

Optional Accessories