

Risk Management Panel PC, models xxxZEUS-198xxxxxxxxxx and xxxZEUS-248xxxxxxxxxx (where x can be 0-9, A-Z, a-z, "-" or blank)

Company:	ONYX HEALTHCARE INC.
Creation Date:	2017-08-09
Abbreviation:	ZEUS-198/248
Project Number:	170503601_3451
Project Description:	The equipment, models xxxZEUS-198xxxxxxxxxx and xxxZEUS-248xxxxxxxxxx (where x can be 0-9, A-Z, a-z, "-" or blank) are Panel PC that are consisted with touch LCD Panel and all electrical components were soldered on PWB and housed in plastic enclosure and powered by an external AC/DC adaptor. The monitor is intended to use for displaying alpha, numerical and graphical data for the medical purpose.
Version Number:	1.00
Version Description:	Workable version for evaluation purpose
Life Cycle Phase:	Design and Development

Table of Content

The Following Documents have been Compiled as a Report:	Page
Definitions and Abbreviations	3
Project Data	6
Version Data	7
Project Team	13
Definition of the Risk Graph	21
Release Criteria	23
Grey Box	26
Risk Analysis	27
Risk Evaluation Before Measures	37
Risk Graph Before Measures	48
Decision about Measures	49
Risk Control	66
Release Protocol of Measures	73
Evaluation After Measures	151
Risk Graph After Measures	181
Risk Graph Before and After Measures	182
Evaluation of Residual Risk	183
Post Production	184
Overview of the Results of the Risk Management	185

	Implemented	Tested	Released
Name	Andrew Wang	Ryan Chen	Mac Lin
Date, Signature	<i>Andrew Wang</i>	<i>Ryan Chen</i>	<i>Mac Lin</i>



Definitions and Abbreviations

Project: 170503601_3451, ZEUS-198/248 | Version: 1.00

Definitions

Word	Explanation
Abnormal Use	Intended act or intended omission of an act by the user or operator of equipment as a result of conduct that is beyond any reasonable means of risk mitigation by the manufacturer.
Adverse Event	Any device-related event that has or may have caused or contributed to a death, or serious injury/deterioration (of patient, user or other person), or has allegedly involved a device malfunction that would be likely to cause or contribute to a death or serious injury (deterioration) if the malfunction were to recur. Such events require notification to the appropriate regulatory agency per defined regulations.
Approval	Agreement that an item or set of actions is complete, has been reviewed, and that the approver concurs or agrees with the outcome and endorses any follow-up actions through a valid signature. To give formal or official sanction to an item.
Design History File	This contains or references the documents that constitute the compilation of records describing the design history of a finished device. The Design History File is not necessarily in a single location. The Design History File covers the technical requirements for Quality System Regulations and the Medical Device Directive.
Direct Harm	Caused by the energy or substance delivered, or not delivered, by/through the product. (i.e. Ventilator)
Harm	Physical injury or damage to the health of people or damage to property or the environment. NOTE: Include unreasonable psychological stress or unwanted pregnancy as part of "damage to the health of people."
Hazard	Potential source of harm.
Hazardous Situation	The possibility of any degree of injury, actual or potential, associated with an event.
Indirect Harm	An injury caused when a product leads an individual to either perform an action or fail to perform an action based on its performance.
Intended Use	Use of a product, process, or service in accordance with the specifications, instructions, and information provided by the manufacturer.
Labeling	Written, printed, electronic, or graphical information included on the device, packaged with the device, accompanied with the device or marketing material, including web-based promotional material. Labeling includes instructions for use, product labels, data sheets, marketing materials/information and promotional materials, installation and service manuals.
Medical Device	used alone or in combination, including the software necessary for its proper application intended by the manufacturer to be used for human beings for the purpose of: <ul style="list-style-type: none"> - Diagnosis, prevention, monitoring, treatment or alleviation of disease, - Diagnosis, monitoring, treatment, alleviation of or compensation for an injury or handicap, - Investigation, replacement or modification of the anatomy or of a physiological process, - Control of conception, and which does not achieve its principal intended action in or on the human body by pharmacological, immunological or metabolic means, but which may be assisted in its function by such means NOTE: The term "medical device" is defined by national law; therefore, consult the laws of the host country.
Medical Intervention	Intervention by a medical professional to help, treat or cure a condition, this includes medical or surgical.

