



IEC 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance	
Report Reference No.	150401702
Date of issue	2015-05-29
Total number of pages	73
CB Testing Laboratory	Victronic Technology Corporation.
Address	4th Fl 130 Ln 235 Baoqiao Rd Xindian Dist New Taipei, 231 Taiwan.
Applicant's name	Onyx Healthcare Inc.
Address	2F., No.135, Ln. 235, Baoqiao Rd., Xindian Dist., New Taipei City 231, Taiwan
Test specification:	
Standard	IEC 60601-1: 2005 + CORR. 1 (2006) + CORR. 2 (2007) EN 60601-1 : 2006/ A11:2011/ A12:2014
Test procedure	CB Scheme
Non-standard test method	None
Test Report Form No.	60601-1TRF_A
Test Report Form Originator	Victronic Technology Corporation
Master TRF	2011/11/29
Copyright © 2010 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved. This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context. If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo shall be removed This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.	
Test item description	Mobile Cart Computer
Trade Mark	 or ONYX Healthcare Inc.
Manufacturer	Onyx Healthcare Inc. 2F., No.135, Ln. 235, Baoqiao Rd., Xindian Dist., New Taipei City 231, Taiwan
Model/Type reference	VENUS-192xxxxxxxxxx and VENUS-222xxxxxxxxxx (Where x can be 0-9, A-Z, "-" or blank)
Ratings	Switching Power Adapter: (Edac Power Electronics Co Ltd, EM11701F) Input: 100-240Vac, 50-60Hz, 2.5-1.0A Mobile Cart Computer: Input rating: 24Vdc, 6.25A



**Battery pack (two in series):
For each rating: 10.8Vdc, Max. 5800mAh**



**DC output (optional): 9Vdc, 1A; 12Vdc, 1A; 24Vdc, 0.5A (can
only use one voltage at the same time)**





VICTRONIC TECHNOLOGY CORPORATION

Testing procedure and testing location:

**Testing location/ address : Victronic Technology Corporation / 4th Fl 130 Ln 235
Baoqiao Rd Xindian Dist New Taipei, 231 Taiwan.**

Tested by (name + signature).. : Ella Chang

Approved by (+ signature) : Daniel Chen


