

Digital Diagnostic Scope

DSC-100

DSC-100, 2MP resolution, Standard



Features

- 3.5" Full Color TFT-LCD
- USB & AV output Full HD live stream
- Desktop charger
- Manual Focus wheel
- Switchable lens for difference use

Specifications

Control unit

Focus	Manual focus
Picture Resolution	1920 x 1080 pixels
LCD Monitor	3.5" TFT LCD
Image Format	JPEG (Photograph) and H.264 (Video)
Interface	Mini USB 2.0, Micro HDMI (1.4), Wi-Fi (2.4G), Bluetooth (4.2)
File Transfer	Mini USB Port, AV out put
Dynamic Video Output	Composite AV out, or USB live video enable from USB port
File Storage	Micro SD Card. Supports 2G to 32GB by FAT32 Format.
Power Source	nput Spec. : 100~240 VAC, 50/60 Hz, 0.2A; Output Spec. : 5V DC, 1.2A
Charging station Input Spec.	5V DC, 1.2 A
Operating Time	3 hours at 2.5 Watt condition
Charging Time	5 hours by DSC 100
Expected service life (defined by manufacturer)	5 years from the date of initial operation *Proper maintenance is necessary.