



Medical Display for Medical Device

MEDDP-415

15" Slim Medical Display



OSD Control Panel



USB for Touch

Power Input

VGA

Audio

HDMI

Features

- 15 " 1024 x 768 Resolution Display
- P-cap touch
- Stylish And Elegant Design

Specifications

DISPLAY

Size	15"
Resolution	XGA 1024 x 768
Max. Colors	16.2M
Brightness (Typ.)	300 nits
Viewing Angle (CR=10)(Typ.)	170° (H)/170° (V)
Contrast Ratio (Typ.)	700:1
Touch Screen	P-cap Touch, Protective Glass
Power Requirement	External power adapter, AC IN: 100~240V, DC 12V/3.33A

SYSTEM

OS Support	Windows® 10, Windows® 11
------------	--------------------------

I/O

Input Signal	HDMI 1.4 x1, VGA x1
USB Port	USB 2.0 Type A x1 for Touch Screen
Audio	Line-in x1
Speaker	Built in 1W speaker x2

MECHANICAL AND ENVIRONMENTAL

Operating Temperature	0°C~40°C (32°F~104°F)
Storage Temperature	-20°C~60°C (-4°F~176°F)
Storage Humidity	10%~90%@35°C, non-condensing
VESA	75mm x 75mm, 100mm x 100mm
Dimension(W/H/D)	352.5 x 276.7 x 45 mm
Packing Size(W/H/D)	402.5x 326.1 x 80.0 mm
Gross Weight	3.7 kg
Net Weight	2.6 kg
Certifications	CE: EN 60601-1-2:2015 + A1:2021 (V4.1), EN 60601-1:2007/A1:2013/A2:2021 (V3.2) FCC: Part 15 Subpart B / Part 18 Subpart C CCC

Ordering Information

- MEDDP-415HNN-A1-1020
15" LCD Monitor.300 nits.VGAX1.HDMIx1.DC 12V Adapter
- MEDDP-415HPN-A1-1020
15" LCD Monitor.300 nits.VGAX1.HDMIx1.DC 12V Adapter.P-CAP touch
- MEDDP-415HPN-A2-1020
15" LCD Monitor.300 nits.VGAX1.HDMIx1.DC 12V Adapter.P-CAP.with CCC certification