


Medica Display for Medical Device

MEDDP-615

15.6" Slim Medical Display



Features

- 15.6" 1920 x 1080 Resolution Display
- Projected Capacitive Multi-touch
- VGA input, HDMI input, and USB Type-C input
- Stylish and elegant design

Specifications

DISPLAY

Size	15.6"
Resolution	FHD 1920 x 1080
Max. Colors	262K
Brightness (Typ.)	400 nits
Viewing Angle (CR=10)(Typ.)	170° (H)/170° (V)
Contrast Ratio (Typ.)	800:1
Touch Screen	Projective capacitive touch screen
Power Requirement	External power adapter, AC IN: 100~240V, DC 12V/3.33A

SYSTEM

OS Support	Windows® 10, Windows® 11
------------	--------------------------

I/O

Input Signal	HDMI 1.4 x1, VGA x1, USB Type C x1
USB Port	USB 2.0 Type B x1 for Touch Screen
Audio	Line-in x1
Speaker	Built in 1W speaker x2

MECHANICAL AND ENVIRONMENTAL

Operating Temperature	0°C~50°C (32°F~122°F)
Storage Temperature	-10°C ~ 60°C (14°F ~ 140°F)
Storage Humidity	10%~90%@35°C, non-condensing
VESA	75mm x 75mm, 100mm x 100mm
Dimension (W/H/D)	387.8 x 232.9 x 38.4 mm
Packing Size (W/H/D)	425.3 x 269.2 x 45.6 mm
Gross Weight	5.20kg (11.46lb)
Net Weight	1.50kg (3.31lb)

Certifications

CE: EN 60601-1-2:2015 + A1:2021 (V4.1),
EN 60601-1:2007/A1:2013/A2:2021 (V3.2)
FCC: Part 18 Subpart C
UL: ANSI/AAMI ES60601-1: 2005 & A1:2012 & A2:2021 (V3.1)
cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1)

Ordering Information

- MEDDP-615HPN-A1-1020
15.6" LCD Monitor, 400 nits, P-CAP touch, VGA x1, HDMI x1, DC 12V Adapter