


## Medical Display for Medical Device

**MEDDP-615**

15.6" Slim Medical Display

**Features**

- 15.6" 1920 x 1080 Resolution Display
- Projected Capacitive Multi-touch
- VGA input, HDMI input, and USB Type-C input
- Stylish and elegant design

**Specifications****DISPLAY**

|                             |                                                       |
|-----------------------------|-------------------------------------------------------|
| Size                        | 15.6"                                                 |
| Resolution                  | FHD 1920 x 1080                                       |
| Max. Colors                 | 262K                                                  |
| Brightness (Typ.)           | 400 nits                                              |
| Viewing Angle (CR=10)(Typ.) | 170° (H)/170° (V)                                     |
| Contrast Ratio (Typ.)       | 800:1                                                 |
| Touch Screen                | Projective capacitive touch screen                    |
| Power Requirement           | External power adapter, AC IN: 100~240V, DC 12V/3.33A |

**SYSTEM**

|            |                          |
|------------|--------------------------|
| OS Support | Windows® 10, Windows® 11 |
|------------|--------------------------|

**I/O**

|              |                                    |
|--------------|------------------------------------|
| Input Signal | HDMI 1.4 x1, VGA x1, USB Type C x1 |
| USB Port     | USB 2.0 Type B x1 for Touch Screen |
| Audio        | Line-in x1                         |
| Speaker      | Built in 1W speaker x2             |

**MECHANICAL AND ENVIRONMENTAL**

|                       |                                                                                                             |
|-----------------------|-------------------------------------------------------------------------------------------------------------|
| Operating Temperature | 0°C~50°C (32°F~122°F)                                                                                       |
| Storage Temperature   | -10°C ~ 60°C (14°F ~ 140°F)                                                                                 |
| Storage Humidity      | 10%~90%@35°C, non-condensing                                                                                |
| VESA                  | 75mm x 75mm, 100mm x 100mm                                                                                  |
| IP Rating             | Front IP65, Rear IP20                                                                                       |
| Dimension (W/H/D)     | 387.8 x 232.9 x 38.4 mm                                                                                     |
| Packing Size (W/H/D)  | 425.3 x 269.2 x 45.6 mm                                                                                     |
| Gross Weight          | 5.20kg (11.46lb)                                                                                            |
| Net Weight            | 1.50kg (3.31lb)                                                                                             |
| Certifications        | CE: EN 60601-1-2:2015 + A1:2021 (V4.1),<br>EN 60601-1:2007/A1:2013/A2:2021 (V3.2)<br>FCC: Part 18 Subpart C |

**Ordering Information**

- **MEDDP-615HPN-A1-1020**  
15.6" LCD Monitor, 400 nits, P-CAP touch, VGA x1, HDMI x1, DC 12V Adapter